

Office of Technology Transfer Product Development Opportunity

Laser Tweezers and Raman Spectroscopy (LTRS)

This innovative combination of Laser Tweezers and Raman Spectroscopy (LTRS) now makes it possible to use one low-powered diode laser for both trapping and Raman excitation in the study of cellular processes and the diagnosis of cellular disorder. The near infra-red wavelength and power-switching operational design allow real-time Raman measurement without damaging the living cell. East Carolina University seeks to commercialize and develop additional applications of this technology.

Commercial Applications

- detect and identify biological organisms
- aid in bio-terrorism defense
- provide rapid and reliable clinical diagnosis
- monitor and detect water or air quality on-site
- aid research in molecular and cellular biology

Benefits

- low cost, high sensitivity, capable of miniturization
- incorporates simple, user-friendly features
- provides real-time Raman measurement
- allows for the study of biochemical change
- allows for capture and manipulation of biological cells and micron-sized particles in fluid

Key Words


- optical tweezers
- raman spectroscopy
- diode laser

figure 1. Experimental setup. A diode laser (DL) beam passes through an interference filter (IF), a pinhole (PH) and is then introduced into an inverted microscope through a dichroic mirror (DM) to form an optical trap. The Raman scattering light of the trapped particles is collected with the same objective, passes through a holographic notch filter (HNF1) and a pinhole (PH), and is then focused onto the entrance of an imaging spectrograph equipped with a CCD detector. M -reflection mirror; L-lens; BS -beam splitter; Obj -objective lens; EP-eyepiece; VC -video camera; TC -Temperature controller; Lamp -Green-filtered xenon illumination light.



The Technology

A compact LTRS system combines the advantages of NIR Raman spectroscopy and optical tweezers for the characterization of single biological cells with a low-power semi-conductor laser. The laser power-switching technique allows low-power trapping (~ 2.0 mW) and high-power Raman excitation (~ 20 mW) of biological cells. In the laboratory, single red blood cells and yeast cells have been trapped and recorded with the LTRS system, and noticeable differences in the Raman spectra of living and dead yeast cells were demonstrated. Using this system, molecular information from in vitro single biological cells can be obtained and used to understand the fundamental cell processes and to diagnose cellular disorder.

Figure 1: Experimental Set-Up



About the Inventors



Dr. Yong-Qing Li
Department of Physics
East Carolina University
Greenville, NC 27858
Phone: 252.328.1858
e-Mail: liy@mail.ecu.edu



Education:

Ph.D. in Physics, Chinese Academy of Sciences, Shanghai, 1989

Representative Publications:

1. C.A. Xie and Y.Q. Li, "Confocal micro-Raman spectroscopy of single biological cells using optical trapping and shifted excitation difference techniques", *Journal of Applied Physics*, 93, 2982-2986 (2003).
2. C.A. Xie and Y.Q. Li, "Raman spectra and optical trapping of highly refractive and nontransparent particles", *Applied Physics Letters*, 81, 951-953 (2002).
3. Y.Q. Li, D. Guzun, and M. Xiao, "Sub-shot noise optical heterodyne detection using amplitude-squeezed light as local oscillator", *Physical Review Letters*, 82, 5225-5228 (1999).
4. Y.Q. Li, P. Lynam, M. Xiao, and P.J. Edwards, "Sub-shot noise laser Doppler anemometry with amplitude-squeezed light", *Physical Review Letters*, 78, 3105-3108 (1997).

Dr. Mumtaz A. Dinno
Department of Physics
East Carolina University
Greenville, NC 27858
Phone: 252.328.6934
e-Mail: dinnom@mail.ecu.edu



Education:

PhD., University of Louisville, 1973

Representative Publication:

1. C.A. Xie, M. A. Dinno, and Y.Q. Li, "Near-infrared Raman spectroscopy of single optically trapped biological cells", *Optics Letters*, 27, 249-251 (2002).

Commercialization Contact

Mark Foley
Office of Technology Transfer
2200 South Charles Boulevard
East Carolina University
Greenville, NC 27858
Phone: 252.328.9549
Fax: 252.328.0799
E-Mail: foleym@mail.ecu.edu



LASER TWEEZERS &
RAMAN SPECTROSCOPY
(LTRS)